

Application News

No. AD-0205

Food Analysis / LCMS-8060

LC-MS/MS Method for Quantitation of Eight Vitamin B in Infant Milk Powder

Yin Ling Chew¹, Eleanor Wai Yi Hor*, Qi Zong Lee**, Noah Luzheng Ong**, Muhammad Averus²; Suwiton², Martin Uli Tua², Zhaoqi Zhan¹ and Jie Xing¹

¹Shimadzu (Asia Pacific) Pte Ltd, Singapore, ²PT Saraswanti Info Genotech, Bogor City, Indonesia

□ Introduction

Water-soluble B vitamins are essential nutrients needed by human body to function properly. Since B vitamins are not produced in adequate amounts in human body, it is important to obtain them from food sources. Low levels of B vitamins in breastfeeding mothers may cause their babies slow growth. Therefore, infant formulas must include proper amounts of nutrients including B vitamins. The Dietary Reference Intake (DRI) and the Food Additive have issued the reference value of infants' daily requirement of B vitamins in infant formula. HPLC and LC/MS/MS methods have used in analysis of B vitamins and additives, either limited to few compounds in HPLC method or required time-consuming sample clean-up method due to the complicated matrix of infant formula. Here, we present a LC-MS/MS method for simultaneous determination of 8 water-soluble B vitamins in infant formula with 5 isotope labelled internal standards, with a simple sample extraction procedure.

□ Experimental

Standard Preparation

Eight water-soluble vitamin standards namely thiamine (B1), riboflavin (B2), nicotinic acid (B3), nicotinamide (B3), pantothenic acid (B5), pyridoxine (B6), biotin (B7) and folic acid (B9) were obtained from Sigma Aldrich. Five isotope labelled internal standards, ¹³C₄-thiamine (IS B1), ¹³C₄, ¹⁵N₂-riboflavin (IS B2), ²H₄-nicotinamide (IS B3), ¹³C₆, ¹⁵N₂-pantothenic acid (IS B5) and ²H₂-biotin (IS B7) were purchased from IsoSciences. Vitamins B1, B3, B5 and B6 were diluted using water while vitamins B2, B7 and B9 were diluted using 0.01M NaOH to form individual stock solutions. All the vitamin standards and internal standards were diluted into working solutions to construct calibration curves.

Sample Preparation

1 gram of infant formula was weighed into a 50 mL polypropylene tube. 10 mL of water was added and vortexed. 30 mL of 1% acetic acid was then added and

vortexed, and centrifuged. The supernatant was filtered using a 0.22 nylon filter and diluted 5 times further using water.

Analytical conditions

Shimadzu LCMS-8060 with heated ESI, and a Shim-pack GIST C18-AQ (2.1 mm x 50 mm; 1.9 µm) column was used in this analysis.

Table 1: Analytical conditions of 8 water-soluble vitamins analysis with LCMS-8060

Column	Shim-pack GIST C18-AQ 1.9µm, 2.1 mm I.D. x 50 mm L
Flow	0.3 mL/min
Mobile phase	A : 5mM Ammonium Formate in MilliQ water 0.1% formic acid B : Methanol
Oven Temp	40 °C
Injection vol	2 µL
Elution mode	Gradient Elution; B% : 0% (0.0 to 2.0 min) → 30% (2.5 min) → 50% (5.0 min) → 99% (5.5 to 7.0 min) → 0% (7.1 to 12.0 min)
Interface	ESI
MS Mode	Positive mode
CID gas	Argon, 230 kPa
Heat Block Temperature	400 °C
DL Temperature	250 °C
Interface Temperature	300 °C
Nebulizing Gas Flow	Nitrogen, 3.0 L/min
Drying Gas Flow	Nitrogen, 10.0 L/min
Heating Gas Flow	Zero Air, 10.0 L/min

MRM method for eight water-soluble B vitamins and internal standards

Automated MRM optimization of 8 water-soluble vitamins standards and 5 internal standards were

*Student for internship training program from *Nanyang Technological University and **National University of Singapore

Table 2: MRM transitions of 8 water-soluble vitamins (B) and 5 internal standards on LCMS-8060

Vitamin (Compound Name)	Quantifier ion (m/z)	Qualifier ion (m/z)	Internal standard	Quantifier ion	Qualifier ion
B1 (Thiamine)	265.00>81.10	265.00>144.10	IS B1	269.00>122.10	269.00>148.10
B2 (Riboflavin)	377.00>243.15	377.00>172.20	IS B2	383.00>249.15	383.00>175.15
B3 (Nicotinic acid)	124.00>53.10	124.00>78.10	IS B3	127.00>84.10	127.00>81.10
B3 (Nicotinamide)	123.00>53.10	123.00>78.05			
B6	170.00>77.10	170.00>152.10	IS B5	224.00>206.10	224.00>188.10
B5	220.10>90.10	220.10>202.25			
B7	245.00>227.10	245.00>97.10	IS B7	247.00>229.10	247.00>99.10
B9	442.10>295.10	442.10>176.10			




Figure 1: Total ion chromatograms (TIC) of water-soluble B vitamin standards at 50 ng/mL.

carried out using the LabSolutions workstation. The precursors were their protonated ions, [M+H]⁺. Two MRM transitions for each compound were chosen as quantifier and qualifier ions (Table 2). Linear calibration curves were obtained for all 8 water-soluble vitamin compounds. Good linearity with correlation coefficient (r²) greater than 0.999 across the range of 10 ng/mL – 750 ng/mL was obtained. The calibration curves of vitamin B2, B3 and B5 are shown in Figure 2.

Performance evaluation of established quantitation method

The limit of detection (LOD) and limit of quantitation were obtained from LabSolutions software with S/N~3 and S/N~10. The LOD and of LOQ of 8 water-soluble B vitamins ranges from 0.08 – 1.36 ng/mL and 0.23 – 4.13 ng/mL, respectively. Repeatability (n=6) tests were performed at concentration of 10 ng/mL for each vitamin standard and the results are tabulated in Table 3.




Figure 2: Selected calibration curves of water-soluble vitamin standards with internal standards.

Table 3: Linearity, LOD, LOQ and repeatability results for 8 water-soluble vitamin compounds

Vitamin	Range (ng/mL)	R2	LOQ (ng/mL)	LOD (ng/mL)	RSD% (n=6)
B1	1-500	0.9994	0.52	0.17	2.80
B2	1-750	0.9999	0.23	0.08	3.67
B3 (Nicotinic acid)	1-750	0.9993	0.97	0.32	3.20
B3 (Nicotinamide)	1-750	0.9994	0.52	0.17	2.36
B5	5-500	0.9990	4.13	1.36	1.91
B6	1-750	0.9999	0.51	0.17	3.94
B7	1-100	0.9991	0.32	0.11	2.07
B9	1-100	0.9992	1.00	0.33	2.87

Recovery study of 8 water-soluble vitamins in infant formula powder

A recovery test was performed using infant formula powder as the matrix. The 8 water-soluble vitamins are spiked in the infant formula powder at different concentrations namely, 25 mg/kg for vitamins B3 (both nicotinic acid and nicotinamide) and B5; 5 mg/kg for vitamins B1, B2 and B6; 0.5 mg/kg for vitamins B7 and B9.

The recovery test was performed twice to obtain reliable results. Each duplicate was injected thrice and the average concentration was used for the calculation of recovery. Good recovery of 8 water-soluble vitamins was achieved between 90.1~ 113.6%.

Quantitation of 8 water-soluble vitamins in NIST 1849a

The established method was evaluated using the standard reference material, NIST 1849a, which is an infant/adult nutritional formula. Two sets of NIST 1849a milk powder was extracted and each of them was injected in triplicates to ensure the accuracy of the method. The results obtained are within the range of content and are tabulated in Table 4.

Table 4: Average recovery of water-soluble vitamins in infant formula (n=2 x 3).

Compound name	Recovery (%)
B1	113.6
B2	105.2
B3 (Nicotinic Acid)	101.6
B3 (Nicotinamide)	109.8
B5	104.9
B6	108.0
B7	90.1
B9	90.1

Table 5: Results of quantitation of 8 B vitamins in NIST 1849 nutritional formula (D.F. = 200).

Vitamin (Compd. name)	Content (mg/kg)	Replicate (mg/kg) (n=3)		
		1st	2nd	3rd
B1 (Thiamine)	12.57 ± 0.98	13.14	13.23	12.89
B2 (Riboflavin)	20.37 ± 0.52	20.59	20.80	20.77
B3 (Nicotinamide)	108 ± 10	100.92	101.31	100.53
B5 (Pantothenic acid)	68.2 ± 1.9	69.02	67.99	67.34
B6 (Pyridoxine)	13.46 ± 0.93	13.58	14.03	13.82
B7 (Biotin)	1.99 ± 0.13	2.06	2.09	2.09
B9 (Folic acid)	2.293 ± 0.062	2.253	2.285	2.277




Figure 3: Selected MRM chromatograms of NIST 1849a sample

Conclusions

A sensitive and selective LC-MS/MS method with a simple extraction procedure was established for quantitation of 8 water-soluble B vitamins using 5 isotope labelled internal standards. The method performance of analysis and sample preparation were evaluated. The linearity (r^2) of all the 8 B vitamins are greater than 0.999. The LOD and LOQ estimated are at 0.08~1.36 ng/mL and 0.23~4.13 ng/mL, respectively. The recovery ranges from 90.1% to 113.6%. The method was used to quantitate the standard reference material, NIST 1849a. The quantitation results fall within the official values given.

References

- Cellar, N. A.; McClure, S. C.; Salvati, L. M.; Reddy, T. M. *Anal. Chim. Acta* . **2016**, 934, 180–185.
- Goldschmidt, R. J.; Wolf, W. R. *Anal. Bioanal. Chem.* **2010**, 397 (2), 471–481.